# Series 2400 Total Hydrocarbon Analyzer

GOW-MAC

Operates on 100%  $\rm H_{_2},$  40/60%  $\rm H_{_2}/\rm N_{_2}$  or H\_/He Fuel Mixes

Methanizer to measure CO and  $CO_2$ 

Consistent, Precise, Flexible, Compact

19" Rack Mount or Bench Top

The GOW-MAC Series 2400 microprocessor-controlled total hydrocarbon analyzer is designed to continuously measure concentrations of hydrocarbons and, utilizing an optional catalytic methanizer, CO and  $\rm CO_2$  in gas streams.

## Principal of Operation

The flame ionization detector (FID) operates on the principle that the electrical conductivity of a gas is directly proportional to the concentration of charged particles within the gas. Ionized carbon atoms are produced when burned in a hydrogen flame. Positive ions, negative ions and electrons present between an electrode gap decrease the gap resistance and cause an electrical current to flow across the gap.<sup>1</sup> The resulting voltage decrease is detected and then reported as a concentration.

### Applications

The instrument provides long term stability and accurate performance in a multitude of gas manufacturing facilities and laboratory applications, such as:

- · Gas purity certification
- · UHP inert gas contaminant analysis
- · Carbon bed breakthrough detection
- HC monitoring in CO<sub>2</sub> feed streams in the food/ beverage industry
- Process monitoring
- · Monitoring ambient air and other gases
- Monitoring of HC emissions
- R&D laboratory applications

## Features

- Flame Ionization Detector (FID)
- Methanizer to measure CO and CO<sub>2</sub>
- High resolution graphical LCD touchscreen
- · Flameout indicator / Auto shutoff for gases
- Auto ranging from 0.00 ppm to 10,000 ppm
- · Auto zero and calibration

Series 2400 Total Hydrocarbon Analyze

GOW-MAC

- · Programmable relays for concentration alarms
- Electronic flow control of air, fuel, and sample gas
- · Real-time data logging by USB and serial port
- Outputs: Analog (standard): 4-20 mA, VDC [0-5 V] or Digital (optional): RS-232

The Series 2400 analyzer incorporates a flow system that can be used for ambient air analysis where low levels of hydrocarbons are monitored on a continuous basis. Samples can be introduced by pressurized tanks or gas lines.

#### Flow Control

Electronic pressure controllers (EPCs) are used to deliver regulated streams of air, fuel and sample to the FID.

<sup>1</sup> McNair, H.M.; Bonelli, E.J. Basic Gas Chromatography; Varian; Berkeley, CA 1968

## Specifications

Detector:	Flame Ionization (FID)	
Sensitivity:	$0.05 \text{ ppm (CH}_4)$	
Fuel Gas Requirements:	Zero Grade H <sub>2</sub> : 5-100 cc/min clean, zero grade H <sub>2</sub> (THC $\leq$ 0.2 ppm) @ 20 psig Fuel Mix: 5-100 cc/min 40/60% fuel mix (THC $\leq$ 0.5 ppm) @ 20 psig*	
Sample Gas Require- ments:	Minimum of 2 psig with a minimum flow rate of 5-70 cc/min	
Flame Air Requirements:	Minimum of 25 psig with a minimum flow rate of 10-350 cc/min zero air (THC $\leq$ 1.0 ppm)	
Operating Temperature:	mperature: 15 °C to 35 °C (59 °F to 95 °F) [best performance @ stable temperature]	
User Interface 4" W x 2.4 H" graphical resistive touchscreen LCD display (4.7" diagonal)		
Flow Control:	Electronic pressure control (EPC) of fuel, air, and sample	
Operation:	Automatic zero at user-defined intervals	
Ranges:	4 Ranges: (0-10, 0-100, 0-1,000, 0-10,000) or Autorange	
Inputs:	Span (calibration)	
Outputs:	Analog: 4-20 mA, VDC [0-5 V] (standard) Digital: USB and Serial (Data Logging)	
Relays:	Programmable relays for two (2) high/low concentration alarms [2 1A alarm relays and 1 ready relay]	
Safety:	Flame-out Circuit with automatic shut-off of fuel	
Response Time:	<5 seconds to 90% of final reading	
Drift:	< 1% full scale drift over 24 hours	
Reproducibility:	± 1% full scale for successive samples	
Mounting:	3U case, 19" rack mount or bench top	
Connections:	1/8" tube fitting connectors	
Power:	Series 2400, -CM: 100 W @ 115 V, 60 Hz, 1.0 amps Series 2402, -CM: 100 W @ 230 V, 50 Hz, 0.5 amps	
Dimensions:	19.00" W x 5.25" H x 20.50" D (48.3 cm x 13.3 cm x 52.1 cm)	
Weight (typical):	Net: 35lbs. (15.9 kg)	

\*Note: There is a decrease in sensitivity when using a fuel mix other than zero grade H<sub>2</sub>. Fuel mix required for oxygen matrix.

Specifications and features will vary depending upon system configuration and are subject to change without notice. The above specifications are established during design, but are not to be construed as test criteria for every product.

## **Ordering Information**

$115 \mathrm{V}$	$230 \mathrm{V}$	
<u>Model</u>	<u>Model</u>	
2400	2402	Total Hydrocarbon Gas Analyzer, standard
2400-CM	2402-CM	Total Hydrocarbon Gas Analyzer with catalytic methanizer

Global Headquarters

277 Brodhead Road

Tel: (610) 954-9000

Fax: (610) 954-0599

GOW-MAC Instrument Co.

Bethlehem, PA 18017 U.S.A.

E-mail: sales@gow-mac.com URL: www.gow-mac.com





SB-2400 0921

©2014GOW-MAC Instrument Co. All rights reserved. The GOW-MAC logo and design are registered trademarks of GOW-MAC Instrument Co. in the United States and other countries. Information, descriptions, photographs and specifications in this publication are subject to change without notice. GOW-MACInstrument Co. shall not be liable for errors contained here in or for incidental or consequential damages in connection with the furnishing, performance or use of this material.

Asian-Pacific Office

Tel: +886-3-5422845

Fax: +886-3-5422045 E-mail: jclo@gow-mac.com

GOW-MAC Instrument Co. - Taiwan 3Fl. No. 33, Lane 92, Beida Road

Hsinchu City 300, Taiwan ROC

China Office

Beijing, China

Tel: + 86-13-501010846

GOW-MAC Instrument Co. - China

E-mail: sales@gow-mac.com.cn

URL: www.gow-mac.com.cn

